

H160ZFP Foot Pedal System

Remote Focus Control Units For Most Optical Systems

Features



Prior Scientific's H160ZFP Foot Pedal System provides a cost-effective remote focusing solution for a variety of applications.

Features of the H160ZFP:

- Units can be fitted on virtually any new or existing microscope with coaxial coarse and fine focus
- Improves ergonomics
- Reduces the risk of repetitive stress injuries
- Units are simple and easy to operate
- Provides quick, precise, maintainable focus at various speeds
- Small footprint allows optimal placement on desk or bench top
- Economical design provides a cost-effective remote focusing solution
- Provides hands free operation for surgical procedures

All H160ZFP Foot Pedal Systems include a control box with on/off switch, speed control, drive motor, foot pedals, focus mount and focus sleeve to match your new or existing microscope.

Prior products have a well deserved reputation for quality and repeatability. For the remote focus control units and all of their products, Prior provides full support and service - both direct and indirect - through a professional, knowledgeable and extensive dealer network.

Ordering Information

Part Number: H160ZFP	Foot pedal control DC drive remote focus control unit for microscopes with coaxial coarse and fine focus knobs
-------------------------	--

PRIOR SCIENTIFIC INC.,
80 RESERVOIR PARK DRIVE
ROCKLAND, MA 02370-1062
TEL 781-878-8442 FAX 781-878-8736
WWW.PRIOR.COM