


Micromanipulator Mounting Plate

For BXWI and IX Series Flat Top Stages

Features



PRIOR
Scientific

Prior Scientific is proud to introduce the HPI 139 Micromanipulator Mounting Plate for both the BXWI and the IX series flat top stages. Providing access from the front, left and right sides of the stage, the HPI 139 mounting plate maximizes the area available for mounting micromanipulators. The HPI 139 comes with the standard 1/4-20 UNC tapped holes on 1" centers.

Ordering Information

Part #: HPI 139	Micromanipulator Mounting Plate for BXWI & IX Series Flat Top Stages
-----------------	--