

## PL220 Automated Oiler

An automated system to add immersion oil to microscope slides



The PL-220 Automatic Oiler takes screening to the next level by allowing for the use of high resolution, high numerical aperture oil immersion objectives within an automated system.

- Allows even greater levels of automation and flexibility within optical microscopy.
- Provides highly controlled automatic delivery of immersion oil onto the slide surface.
- Easy to use
- Simple software integration
- Compatible with all leading upright microscopes.
- Stand alone control unit
- USB (HID or Virtual COM) or with RS232C interfaces
- Suitable for customised or OEM applications

### Products available

Product	Description
PL220	Automatic Oil Dispenser. A programmable system to automatically add immersion oil to microscope slides. Includes mounting brackets and reservoir.

### Worldwide distribution

**Prior Scientific Ltd**  
Cambridge, UK  
T. +44 (0) 1223 881 711  
E. [uksales@prior.com](mailto:uksales@prior.com)

**Prior Scientific Inc**  
Rockland, MA USA  
T. +1 781-878-8442  
E. [info@prior.com](mailto:info@prior.com)

**Prior Scientific GmbH**  
Jena, Germany  
T. +49 (0) 3641 675 650  
E. [jena@prior.com](mailto:jena@prior.com)

**Prior Scientific KK**  
Tokyo, Japan  
T. +81-3-5652-8831  
E. [info-japan@prior.com](mailto:info-japan@prior.com)