

HF108IX3 Filter Wheel

For Olympus IX3 Inverted Microscope Systems



Prior Scientific is pleased to introduce the new HF108IX3 Filter Wheel deck designed for the Olympus IX3 inverted microscope systems. Simple cassette like insertion into the side deck of the Olympus IX3 makes installation quick and easy.

- The HF108IX3 provides a compact high speed filter wheel option which accommodates up to eight 32 mm diameter filters
- 32mm filters completely fill the field of view in the eyepieces and large chip cameras
- Filter wheel mounted below the nosepiece allows the user to not only switch between wavelengths at high speed, but also to see the fluorescent images through the eyepieces of the IX3 microscope as well as the camera
- Compatible with ZDC (Z Drift Compensator)
- ProScan and OptiScan integration as standard filter wheel
- Software compatible with all Prior products
- Drivers available for cellSens and Metamorph

General Specifications:

- · Capacity of eight 32mm diameter filters
- Compatible with ProScan and OptiScan controllers
- · Drivers available for cellSens and Metamorph
- Compatible with ZDC (Z Drift Compensator)



HF108IX3 Filter Wheel

Complete System including H31ADF Module (requires existing ProScan[®]III H31XYZE Controller)

Part Number	Description
HF108IX3B	HF108IX3 filter wheel and H31ADF control module

IX3 Filter Wheel

Part Number	Description
HF108IX3	High speed filter wheel deck module for Olympus IX3 microscope

IX3 Filter Wheel Control Options

Part Number	Description
ProScan	Compatible with all ProScan controllers which include a filter wheel drive
OptiScan	Compatible with all OptiScan controllers which include a filter wheel drive

Sample Fluorescent System Configuration for Typical Olympus IX83 Microscope

Part Number	Description
IX3PROXY	ProScan XY System for IX83, Including ProScan stage, ProScan III controller, joystick and universal sample holder
HF108IX3B	HF108IX3 filter wheel and H31ADF control module
HF110	Excitation filter wheel
HF235	Excitation adapter IX3
L200SIX2	Lumen 200S fluorescence illuminator with internal high speed shutter
LDB100OL	Brightfield LED for Olympus models





Worldwide distribution

Prior Scientific, Ltd Cambridge, UK

T. +44 (0) 1223 881711 E. uksales@prior.com

Prior Scientific, Inc Rockland, MA. USA

T. +1 781-878-8442 E. info@prior.com

Prior Scientific, GmbH Jena, Germany

T. +49 (0) 3641 675 650 E. jena@prior.com

Prior Scientific KK Tokyo, Japan

T. +81-3-5652-8831 E. info-japan@prior.com