



## ES10ZE Controller

Motor controller for focus control



This Focus only controller from Prior offers an affordable solution for microscope automated focus. Compatible with a wide range of microscopes and third party imaging software, the ES10ZE offers fine and precise control over focussing.

With a compact design of only 160 x 125 x 60 mm the ES10ZE takes up little valuable laboratory space, and has both RS232 and USB connections, increasing its versatility. A focus encoder input is also present.

The controller also has a manual focus knob, and the screen offers an accurate position readout. Convenience of control is greatly increased by the presence of a number of speed options, and fast 'Up' and 'Down' options.



**Prior Scientific Ltd**

3-4 Fielding Industrial Estate • Wilbraham Road • Fulbourn • Cambridge • CB21 5ET • UK  
t: +44 (0)1223 881711 • e: [uksales@prior.com](mailto:uksales@prior.com) • [www.prior-scientific.co.uk](http://www.prior-scientific.co.uk)



FM 61600



**Prior Scientific Inc**

80 Reservoir Park Drive • Rockland • MA. 02370 • U.S.A.  
t: +1 781-878-8442 • e: [info@prior.com](mailto:info@prior.com) • [www.prior-us.com](http://www.prior-us.com)

Registered Address: Units 3/4 Fielding Industrial Estate • Wilbraham Road • Fulbourn • Cambridge • CB21 5ET • United Kingdom Registered in England 404087