


## LDBUSBTTL Adapter

For USB to TTL output for simple device control



The LDBUSBTTL adapter is a USB to TTL output converter from Prior Scientific that allows users to control simple devices such as TTL controlled shutters and transmitted LED light sources.

- Powered via USB
- Includes a USB cable and a converter box with a female BNC.
- Plug and play functionality as a Prior controller.
- Maintained output of 0V/5V

The LDBUSBTTL adapter is compatible with:

- Prior LDB100/101F Brightfield LED
- Lumen 200S
- Third party TTL controlled shutters

Note that this device is for TTL output only (from PC to external device).

### Products available

Product	Description
LDBUSBTTL	USB to TTL output converter
W3046	USB Cable - included with LDBUSBTTL

### Worldwide distribution