


W4416/PLW110 4 Slide Holder in Well Plate Format



The W4416 slide holder can hold up to four conventionally sized (25mm x 75mm) microscope slides in a well plate format.

The holder is stackable, minimising the space taken up, and a simple 'lock in' mechanism allows the end user to quickly and easily insert the slides. Once inserted, the slides are held safely and securely in place - the W4416 can even be inverted and the slides will remain in place.

The W4416 slide holder will prove to be useful for researchers using inverted microscopes to image conventional slides, especially when large numbers of samples must be rapidly imaged.

Ordering Information

Code	Description
W4416	4 slide holder in wellplate format (x1)
PLW110	4 slide holder in wellplate format (x10)

Worldwide distribution

Prior Scientific Ltd
Cambridge, UK
T. +44 (0) 1223 881 711
E. uksales@prior.com

Prior Scientific Inc
Rockland, MA USA
T. +1 781-878-8442
E. info@prior.com

Prior Scientific GmbH
Jena, Germany
T. +49 (0) 3641 675 650
E. jena@prior.com

Prior Scientific KK
Tokyo, Japan
T. +81-3-5652-8831
E. info-japan@prior.com