

ES10ZE Controller

Motor controller for focus control



This Focus only controller from Prior offers an affordable solution for microscope automated focus. Compatible with a wide range of microscopes and third party imaging software, the ES10ZE offers fine and precise control over focussing.

With a compact design of only 160 x 125 x 60 mm the ES10ZE takes up little valuable laboratory space, and has both RS232 and USB connections, increasing its versatility. A focus encoder input is also present.

The controller also has a manual focus knob, and the screen offers an accurate position readout. Convenience of control is greatly increased by the presence of a number of speed options, and fast 'Up' and 'Down' options.

Worldwide distribution

Prior Scientific Ltd
Cambridge, UK
T. +44 (0) 1223 881 711
E. uksales@prior.com

Prior Scientific Inc
Rockland, MA USA
T. +1 781-878-8442
E. info@prior.com

Prior Scientific GmbH
Jena, Germany
T. +49 (0) 3641 675 650
E. jena@prior.com

Prior Scientific KK
Tokyo, Japan
T. +81-3-5652-8831
E. info-japan@prior.com