



DIMENSIONS ARE SHOWN FOR H117P1N4.

MOST H117 STAGES FOR NIKON MICROSCOPES MOUNT IN THE IDENTICAL FASHION.

DIMENSIONS ARE IN MILLIMETERS
 TOLERANCE
 ANGULAR: +/- 1/2 Degree
 ONE PLACE DECIMAL ± 0.4MM
 TWO PLACE DECIMAL ± 0.1MM
 REMOVE BURRS AND
 BREAK SHARP EDGES
 DO NOT SCALE DRAWING

NAME	DATE
ENG APPR. JRD	5/5/2023
PROPRIETARY INFORMATION OF:	
PRIOR SCIENTIFIC INC 80 RESERVOIR PARK DRIVE ROCKLAND, MA 02370	

PRIOR SCIENTIFIC INC		
TITLE: H117P1N4 MECHANICAL OUTLINE		
SIZE B	DWG. NO. H117P1N4-OUTLINE	REV A
SCALE: 1:2		SHEET 1 OF 1